

 PRODUCT  
INFORMATION

Liquid Handling Station

## BRAND Liquid Handling Station


Contamination-free, automatic pipetting into 384-well microplates with the BRAND Liquid Handling Station

### Introduction

With the BRAND Liquid Handling Station (LHS), a wide variety of everyday laboratory tasks can be performed automatically, with accurate results. Manually pipetting into 384-well microplates, as is common in high-throughput screening (HTS), is easily susceptible to errors and contamination. The objective of this technical note is to determine the contamination risk for pipetting into a 384-well microplate with the LHS. A common method used to establish this risk is to pipette a measuring solution into a microplate in a checkerboard pattern. Because the measuring solution is missing in every second well, contamination between wells can be measured.

### Material and method

Every second well of a 384-well microplate (Cat. No. 781620), beginning with Well A2, is filled in checkerboard pattern with 50  $\mu$ l of a 1.5% solution of the food coloring Patent Blue V (E131). Next, each empty well is filled with 50  $\mu$ l of fully deionized (DI) water. The pipetting procedure is carried out once with the 50  $\mu$ l single-channel liquid end (SC LE) and once with the 50  $\mu$ l multi-channel liquid end (MC LE). In addition, in another plate, each well is filled with 50  $\mu$ l of DI water. This plate serves as the negative control. The absorption of all plates is then measured in the Nanoquant Infinite M200 Pro from Tecan at 640 nm. Based on the signal strengths in the wells, in which only DI water had been pipetted, it can then be determined if contamination of the food coloring has taken place.



384-well microplate (781620)



**BRAND GMBH + CO KG**

P.O. Box 1155 | 97861 Wertheim | Germany

T +49 9342 808 0 | F +49 9342 808 98000 | [info@brand.de](mailto:info@brand.de) | [www.brand.de](http://www.brand.de)



**BRAND. For lab. For life.®**

BRAND®, BRAND. For lab. For life.®, as well as the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

California Residents: For more information concerning California Proposition 65, please refer to [www.brand.de/calprop65](http://www.brand.de/calprop65).

Subject to technical modification without notice. Errors excepted.

© 2021 BRAND GMBH + CO KG | Printed in Germany | 0621



Find accessories and replacement parts, user manuals, test instructions (SOP) and product videos at [shop.brand.de](http://shop.brand.de)



Further information on products and applications can be found on our YouTube channel: [mylabBRAND](#)

BRAND (Shanghai) Trading Co., Ltd.  
Shanghai, China

Tel.: +86 21 6422 2318  
[info@brand.com.cn](mailto:info@brand.com.cn)  
[www.brand.cn.com](http://www.brand.cn.com)

BRAND Scientific Equipment Pvt. Ltd.  
Mumbai, India

Tel.: +91 22 42957790  
[customersupport@brand.co.in](mailto:customersupport@brand.co.in)  
[www.brand.co.in](http://www.brand.co.in)

BrandTech® Scientific, Inc.  
Essex, CT. United States of America

Tel.: +1 860 767 2562  
[info@brandtech.com](mailto:info@brandtech.com)  
[www.brandtech.com](http://www.brandtech.com)